

# Automatic self zeroing glass buret

Buret for direct reading in °SH based on 25 ml milk sample, automatic self zeroing, Schellbach stripping.

| Type   | PK | Cat. No. |
|--|----|----------|
| Buret: 0 to 25 °SH : 1/10 °SH, polyethylene reservoir 500 ml with base | 1  | 07.0100  |
| Buret: 0 to 40 °SH : 1/5 °SH, polyethylene reservoir 1000 ml with base | 1  | 07.0110  |

- Zubehör:  
Replacement Buret / Gerber Instruments AG  
Test Tube for Reduction Test / Gerber Instruments AG  
Measuring Pipets graduated to Acidity Determination / Gerber Instruments AG  
Test Tube Rack / Gerber Instruments AG