

Test sieves

With welded seam (heat resistant up to 170 °C). Diameter: 200 mm.
Height: 50 mm.

Mesh size µm	Ø mm	Height mm	Material	PK	Cat. No.
20	200	50	Stainless steel	1	9.226 101
25	200	50	Stainless steel	1	9.226 103
32	200	50	Stainless steel	1	9.226 105
36	200	50	Stainless steel	1	9.226 106
38	200	50	Stainless steel	1	9.226 107
40	200	50	Stainless steel	1	9.226 108
45	200	50	Stainless steel	1	9.226 109
50	200	50	Stainless steel	1	9.226 110
53	200	50	Stainless steel	1	9.226 111
56	200	50	Stainless steel	1	9.226 112
63	200	50	Stainless steel	1	9.226 113
71	200	50	Stainless steel	1	9.226 114
75	200	50	Stainless steel	1	9.226 115
80	200	50	Stainless steel	1	9.226 116
90	200	50	Stainless steel	1	9.226 117
100	200	50	Stainless steel	1	9.226 118
106	200	50	Stainless steel	1	9.226 119
112	200	50	Stainless steel	1	9.226 120
125	200	50	Stainless steel	1	9.226 121
140	200	50	Stainless steel	1	9.226 122
150	200	50	Stainless steel	1	9.226 123
160	200	50	Stainless steel	1	9.226 124
180	200	50	Stainless steel	1	9.226 125
200	200	50	Stainless steel	1	9.226 126
212	200	50	Stainless steel	1	9.226 127
224	200	50	Stainless steel	1	9.226 128
250	200	50	Stainless steel	1	9.226 129
280	200	50	Stainless steel	1	9.226 130
300	200	50	Stainless steel	1	9.226 131
315	200	50	Stainless steel	1	9.226 132