

Dry Milk Butyrometer Teichert without Accessories

Weighing method for 2.5 g dry milk, closed bulb with frosted label, neck & Oslash; 17/18 mm

%Fat		PK	Cat. No.
0 - 35	0.25	1	03.07350
0 - 35	0.25	10	03.07350P
0 - 70	0.5	1	03.07700
0 - 70	0.5	10	03.07700P

Zubehör:

Rubber Stoppers for butyrometer with two openings / Gerber Instruments AG

Glass nail for milkpowder butyrometer / Gerber Instruments AG

Weighing Funnel of Glass / Gerber Instruments AG

Alternativen:

Dry Milk Butyrometer Teichert with accessories / Gerber Instruments AG