

Skim Milk Butyrometer

with frosted label on bulb, with or without stopper

%Fat	Division	Milk quantity	Type	PK	Cat. No.
0 - 0.5	ring	43.5 ml (± 0.5)	closed bulb, plain neck, Siegfeld, without stopper	1	03.0211
0 - 0.5	ring	43.5 ml (± 0.5)	closed bulb, plain neck, Siegfeld, without stopper	10	03.0211P
0 - 1	plain	10.73 - 11 ml	open bulb, plain neck, Sichler, with one stopper	1	03.0231
0 - 1	plain	10.73 - 11 ml	open bulb, plain neck, Sichler, with one stopper	10	03.0231P
0 - 1	plain	10.73 - 11 ml	open bulb, plain neck, Sichler, without stopper	1	03.02310
0 - 1	plain	10.73 - 11 ml	open bulb, plain neck, Sichler, without stopper	10	03.02310P

Zubehör:

GERBAL-Stoppers / Gerber Instruments AG

Regulating Pin for GERBAL-Stoppers / Gerber Instruments AG

Rubber Stopper for Milk-, Skim Milk- and Koehler-Cream Butyrometer / Gerber Instruments AG

Rubber Stoppers for butyrometer with two openings / Gerber Instruments AG