

PARAPLAST^{®} Tissue-embedding wax

Granules for embedding tissue. PARAPLASTPLUS^{®} additionally contains dimethyl sulphoxide (DMSO) which improves speed and penetration into tissue samples. Melting point: 55 - 57 °C.

Type	Capacity kg	PK	Cat. No.
PARAPLAST ^{&reg;}	1	1	9.170 015
PARAPLAST ^{&reg;}	5	1	6.255 599
PARAPLASTPLUS ^{&reg;}	1	1	9.170 020