

# Safety pipette bulb Howorka-Ball<sup>&reg;</sup>

- no rigid connection with tube
- no valves or side vent on tubes
- suitable for all standard tube diameters, including pipettes, burettes, filter tubes etc.
- no unstable air pockets
- long operational life due to durable components
- Easily cleaned and steam sterilized if contaminated with liquids

Max. capacity ml	PK	Cat. No.
20	1	9.273 826
50	1	9.273 827

Alternativen:  
LLG-Safety pipette bulb, rubber, red / LLG Labware