

# Automatic self zeroing glass buret

Buret for direct reading in °SH based on 25 ml milk sample, automatic self zeroing, Schellbach stripping.

Type	PK	Cat. No.
Buret: 0 to 25 °SH : 1/10 °SH, polyethylene reservoir 500 ml with base	1	07.0100
Buret: 0 to 40 °SH : 1/5 °SH, polyethylene reservoir 1000 ml with base	1	07.0110

## Zubehör:

Replacement Buret / Gerber Instruments AG

Test Tube for Reduction Test / Gerber Instruments AG

Measuring Pipets graduated to Acidity Determination / Gerber Instruments AG

Test Tube Rack / Gerber Instruments AG