

5-step irreversible temperature strips Thermax^{®}:

5-step irreversible temperature measuring strips covering the temperature range 29 to 290 °C.

Dimensions (WxL): 15 x 35 mm

Type	Min. temperature °C	Max. temperature °C	Temperature	Reversible	Width mm	Length mm	PK	Cat. No.
S	29	40	29 ... 40 °C	No	15	35	10	9.246 310
A	37	46	37 ... 46 °C	No	15	35	10	9.246 311
B	49	71	49 ... 71 °C	No	15	35	10	9.246 312
C	77	99	77 ... 99 °C	No	15	35	10	9.246 313
D	104	127	104 ... 127 °C	No	15	35	10	9.246 314
E	132	154	132 ... 154 °C	No	15	35	10	9.246 315
F	160	182	160 ... 182 °C	No	15	35	10	9.246 316
G	188	210	188 ... 210 °C	No	15	35	10	9.246 317
H	216	249	216 ... 249 °C	No	15	35	10	9.246 318
I	249	290	249 ... 290 °C	No	15	35	10	9.246 319

Alternativen:

8-step irreversible temperature strips testoterm[®] / Testo SE & CO KGaA