

Milk butyrometer Original Gerber

With frosted label on bulb, plain neck, without stopper, for milk quantities of 10.73 - 11 ml.

%Fat	standard	PK	Cat. No.
0 - 4	ISO	1	03.0029
0 - 4	ISO	10	03.0029P
0 - 5	ISO	1	03.0030
0 - 5	ISO	10	03.0030P
0 - 5		1	03.0001
0 - 5		10	03.0001P
0 - 6	ISO	1	03.0031
0 - 6	ISO	1	03.0031P
0 - 7	ISO	1	03.0032
0 - 7	ISO	10	03.0032P
0 - 8	ISO	1	03.0033
0 - 8	ISO	1	03.0033P
0 - 9	ISO	1	03.0009
0 - 9	ISO	10	03.0009P
0 - 10	ISO	1	03.0010
0 - 10	ISO	10	03.0010P
0 - 12		1	03.0012
0 - 12		10	03.0012P
0 - 14		1	03.0014
0 - 14		10	03.0014P
0 - 16		1	03.0016
0 - 16		10	03.0016P
0 - 20		1	03.0020
0 - 20		10	03.0020P
0 - 5 -12		1	03.0025
0 - 5 -12		10	03.0025P

Zubehör:

GERBAL-Stoppers / Gerber Instruments AG

Regulating Pin for GERBAL-Stoppers / Gerber Instruments AG

Rubber Stopper for Milk-, Skim Milk- and Koehler-Cream Butyrometer / Gerber Instruments AG

Alternativen:

Milk butyrometer officially calibrated / Gerber Instruments AG