

# Ultrasonic cleaner baths DU - digital series

With housing made of stainless steel. Timer for countdown (1 to 99 min, model DU-06 3 to 6 min) and built-in heating to 80&deg; C (model DU-06 has no heating). Easy to use with digital controller and simultaneous indication of actual temperature, set temperature, timer and set power. Is possible to use as waterbath with temperature accuracy &plusmn; 3&deg; C. The model DU-100 has an emptying valve.

## Application:

- Cleaning of technical glassware like pipettes, petri dishes and laboratory flasks
- Disinfection and cleaning at the same time
- Degassing of wine or beer samples for analysis of alcohol contents, original worth, color, pH value
- Degassing of food sampled from cans for analysis of stannous contents
- Extraction of soil samples for determination of hydrocarbons

Type	Usable volume (l)	Internal dimension	Outer Dimension	Weight kg	PK	Cat. No.
DU-06	0.6	150 x 95 x 165 mm	160 x 95 x 165 mm	1.35	1	30Gl.702001
DU-32	3.2	240 x 135 x 100 mm	260 x 160 x 235 mm	3.80	1	30Gl.702002
DU-45	4.5	300 x 150 x 100 mm	325 x 175 x 225 mm	4.40	1	30Gl.702003
DU-45	6.5	300 x 150 x 150 mm	325 x 175 x 305 mm	5.40	1	30Gl.702004
DU-100	10	300 x 240 x 150 mm	325 x 265 x 310 mm	8.20	1	30Gl.702005

## Zubehör:

Accessories for ultrasonic cleaning baths AU/DU-series

## Alternativen:

Ultrasonic cleaner baths AU - analogic series