

Centrifuge Gerber Universal

Multi-purpose centrifuge for dairy laboratories.

- 3 pre-programmed speeds for:
- Gerber Butyrometers (350 x g)
 - solubility in milk powder (164 x g)
 - extraction tubes Roese-Gottlieb/Mojonnier (80 x g)
 - 1 variable Speed range (200-1300 rpm, on request)

Microprocessor controller with LED display for speed, run time, temperature and error messages.
With steel housing, locking lid, automatic brake and imbalance detector.
Tempereture adjustable, external dimensions (W x D x H) 640 x 450 x 640 mm.

Rotor	Heating	Voltage	PK	Cat. No.
with Universal-Rotor for 12 positions, without inserts		230 V	1	03.3110-12
with Special-Rotor incl. 36 sleeves		230 V	1	03.3110-36
with Universal-Rotor for 12 positions,without inserts		115 V	1	03.3111-12
with Special-Rotor incl. 36 sleeves		115 V	1	03.3111-36
with Universal-Rotor for 12 positions, without inserts	x	230 V	1	03.3112-12
with Special-Rotor incl. 36 sleeves	x	230 V	1	03.3112-36
with Universal-Rotor for 12 positions, without inserts	x	115 V	1	03.3113-12
with Special-Rotor incl. 36 sleeves	x	115 V	1	03.3113-36

Zubehör:
Rotor for centrifuge Gerber Universal / Gerber Instruments AG
Inserts for Universal rotor / Gerber Instruments AG

Alternativen: