

ASTM Thermometer ACCU-SAFE, calibrated, stem type

Capillary: white back round, environmentally friendly blue filling, with DAkkS-calibration certificate. Immersion depth: total.

Max. Ø: 8 mm
 Colour filling: Blue

Min. temperature °C	Max. temperature °C	Temperature	Division °C	Checkpoints	Length mm	Capillary	Standard
-38	2	-38 ... 2 °C	0.1	3	401	White enameled, round	ASTM S
-20	102	-20 ... 102 °C	0.2	3	435	White enameled, round	ASTM S
-20	102	-20 ... 102 °C	0.2	2	435	White enameled, round	ASTM S
-7	105	-7 ... 105 °C	0.5	3	300	White enameled, round	ASTM S
-7	105	-7 ... 105 °C	0.5	5	300	White enameled, round	ASTM S
-2	80	-2 ... 80 °C	0.2	3	415	White enameled, round	ASTM S
-2	80	-2 ... 80 °C	0.2	1	415	White enameled, round	ASTM S
34	42	34 ... 42 °C	0.1	3	300	White enameled, round	ASTM S
-34	49	-34 ... 49 °C	0.5	3	300	White enameled, prismatic	ASTM S
-8	32	-8 ... 32 °C	0.1	3	401	White enameled, round	ASTM S
-8	32	-8 ... 32 °C	0.1	1	401	White enameled, round	ASTM S