

Hotplate Ceran^{®}

Housing made of stainless steel, brushed and polished, in a unique Softline Design, Schott Ceran^{®} cooking surface. Adjustable power control with power-saving or power control in 9 steps (TC = *Touch Control*). Two-circle HiLight technique, thus extra large heating zone (dia. 210 mm) can be activated. Thermal shutdown and automatic cooking system. Residual heat indicator and warning light. Temperature range (empty): 0 to approx. 460°C. Power supply 230 V.

Type	Dimensions (W x D x H) mm	Hot- plate mm	Heating capacity W	PK	Cat. No.
CT 2200/E	305 x 410 x 75	Ø 140/210	1000/2200	1	9.645 510
CT 2203/TC	305 x 410 x 75	Ø 140/210	1000/2200	1	9.645 511
CT 3400/E	550 x 310 x 65	Ø 140/210, Ø 145	1000/2200, 1200	1	6.261 651