

Ultrasonic cleaner baths AU - analogic series

With housing made of stainless steel. Timer for countdown (1 to 30 min or endless) and built-in heating to 80° C. Analogic control for temperature and timer. AU-450 for degassing of HPLC solvents and processing of samples, with holder SS-200 for sieve cleaning.

Application:

- Cleaning of technical glassware like pipettes, petri dishes and laboratory flasks
- Disinfection and cleaning at the same time
- Degassing of wine or beer samples for analysis of alcohol contents, original worth, color, pH value
- Degassing of food sampled from cans for analysis of stannous contents
- Extraction of soil samples for determination of hydrocarbons

Type	Usable volume (l)	Internal dimension	Outer Dimension	Weight kg	PK	Cat. No.
AU-32	3.2	240 x 135 x 100 mm	260 x 160 x 235 mm	3.80	1	30GI.701001
AU-65	6.5	300 x 150 x 150 mm	325 x 175 x 305 mm	5.40	1	30GI.701002
AU-220	22	500 x 300 x 150 mm	530 x 330 x 340 mm	13.20	1	30GI.701004
AU-450	45	500 x 300 x 300 mm	530 x 325 x 500 mm	17.00	1	30GI.701003

Zubehör:

Accessories for ultrasonic cleaning baths AU/DU-series

Alternativen:

Ultrasonic cleaner baths DU - digital series