

Precision thermometer ACCU-SAFE, similar ASTM, Celsius scale, stem type, immersion

Capillary: colourless prismatic round, diam. 6 to 8 mm, calibratable, environmentally friendly blue filling. With auxiliary scale at 0 °C.

Length: 320 mm
Colour filling: Blue
Capillary: White enameled, round
Standard: ASTM S128C (similar)

Min. temperature °C	Max. temperature °C	Temperature	Division °C	Max. Ø mm	Shape thermometer	PK	Cat. No.
-1,4	1,4	-1,4 ... 1,4 °C	0.05	8	Stem shape	1	4.679 406

Alternativen:

Precision thermometer ACCU-SAFE, similar ASTM, Celsius scale, stem type / Ludwig Schneider

Precision thermometer ACCU-SAFE, similar ASTM, Celsius scale, stem type, immersion depth total / Ludwig Schneider