

Tubing, FEP

FEP fluorocarbon resin. Highly chemical resistant and temperature resistant up to 205°C.

Temperature resistance:

-200 ... 205 °C

Colour:

Transparent

Material:

FEP

Int. Ø mm	Ext. Ø mm	Wall thickness mm	Min. temp. resistance °C	Max. temp. resistance °C	PK	Cat. No.
0.8	1.6	0,40	-200	205	1	9.205 850
1.6	3.2	0,80	-200	205	1	9.205 851
2	4	1,00	-200	205	1	9.205 852
3.6	6	1,20	-200	205	1	9.205 853
4	6.4	1,20	-200	205	1	9.205 854
4	6	1,00	-200	205	1	9.205 855
4.4	6.4	1,00	-200	205	1	9.205 856
5.6	8	1,20	-200	205	1	9.205 857
6	8	1,00	-200	205	1	9.205 858
6.4	9.6	1,59	-200	205	1	9.205 859
6.8	10	1,60	-200	205	1	9.205 860
8	10	1,00	-200	205	1	9.205 861
9.5	12.7	1,59	-200	205	1	9.205 863
10	12	1,00	-200	205	1	9.205 864

Alternativen:

Tubing, PTFE / Bohlender

Tubing, PFA / Bohlender