

Pocket Thermometer

Enclosed scale type, with nickel-plated, double plug collar, ca. 11.5 x 140 mm, with clip, prismatic measuring capillary, uncoated, division: 1 °C.

Type	PK	Cat. No.
measuring range: -40 to +50 °C, red filling	1	10.1205
spare thermometer, measuring range: -40 to 50 °C, red filling	1	10.1207
measuring range: -20 to +120 °C, mercury filling	1	10.1211
measuring range: -20 to +120 °C, blue filling	1	10.1212

Zubehör: