

Sediment Tester

sediment tester Gerber

Foldable aluminum bottle 500 ml, with drain hose and clamp for attachment to the edge of the table, dirt box & Oslash; 18 mm, for filter plates and cotton filters

Type	PK	Cat. No.
Sediment tester Gerber	1	06.8010
Pressure bulb for sediment testers	1	06.8030

Zubehör:

Filters for dirt test in milk / Gerber Instruments AG

Comparison Table for Sediment Test / Gerber Instruments AG

Alternativen:

Sediment Test Pump / Gerber Instruments AG