

Culture flasks, borosilicate glass 3.3, with beaded rim

Erlenmeyer. With 3 side baffles.

Nominal capacity ml	Neck Ø mm	Ø mm	Height mm	Dimensions (ØxH)	No. of baffles	With beaded rim	Material	PK	Cat. No.
250	34	85	140	85 x 140 mm	3	Yes	Boro 3.3	1	9.010 125
300	34	87	156	87 x 156 mm	3	Yes	Boro 3.3	1	9.010 126
500	34	105	175	105 x 175 mm	3	Yes	Boro 3.3	1	9.010 127
1000	42	131	220	131 x 220 mm	3	Yes	Boro 3.3	1	9.010 128