

# Cork borers

Nickel-plated brass. With handles. Tube wall 0.4 mm.

Description	Material	PK	Cat. No.
3-piece set, 5 ... 7.5 mm tube dia.	Brass	1	6.269 657
6-piece set, 5 ... 11.25 mm tube dia.	Brass	1	9.143 106
9-piece set, 5 ... 15 mm tube dia.	Brass	1	9.143 109
12-piece set, 5 ... 18.75 mm tube dia.	Brass	1	9.143 112
18-piece set, 5 ... 26.25 mm tube dia.	Brass	1	9.143 118

Zubehör:

Cork borer sharpener / Carl Friedrich Usbeck KG