

Temperature probe with handle for Handheld thermometer Type 13100

Mineral-insulated thermocouple 1 x NiCr-Ni type K acc. to DIN EN 60584 class 2.
Cable length: 2 m

Min. temperature °C	Max. temperature °C	Temperature	Sort of probe	Type of probe	Ø mm	Length mm	Dimensions (ØxL)
-100	200	-100 ... 200 °C	Immersion probe	Type K	3	150	3 x 150 mm
-100	400	-100 ... 400 °C	Immersion probe	Type K	3	250	3 x 250 mm
-100	800	-100 ... 800 °C	Immersion probe	Type K	3	400	3 x 400 mm
-100	200	-100 ... 200 °C	Insertion probe	Type K	3	150	3 x 150 mm
-100	250	-100 ... 250 °C	Insertion probe	Type K	3	250	3 x 250 mm
-100	250	-100 ... 250 °C	Insertion probe	Type K	3	400	3 x 400 mm
-50	200	-50 ... 200 °C	Insertion probe	Type K	3.8	150	3.8 x 150 mm