

Tubings TYGON^{®} HC F-4040-A, with 2 colour-coded bridges

For Hydro-Carbons and Petro-Chemicals.

- Resists swelling and cracking
- Ozone and UV resistant

Int. Ø mm	Ext. Ø mm	Wall thickness mm	Colour	PK	Cat. No.
0.19	1.99	0.90	Orange/red	12	9.208 302
0.38	2.180	0.90	Orange/green	12	9.208 304
0.51	2.31	0.90	Orange/yellow	12	9.208 305
0.64	2.44	0.90	Orange/white	12	9.208 306
0.76	2.46	0.85	Black/black	12	9.208 307
0.89	2.59	0.85	Orange/orange	12	9.208 308
1.02	2.72	0.85	White/white	12	9.208 309
1.140	2.84	0.85	Red/red	12	9.208 310
1.3	3	0.85	Grey/grey	12	9.208 311
1.42	3.12	0.85	Yellow/yellow	12	9.208 312
1.52	3.22	0.85	Yellow/blue	12	9.208 319
1.65	3.35	0.85	Blue/blue	12	9.208 313
1.85	3.55	0.85	Green/green	12	9.208 314
2.06	3.76	0.85	Purple/purple	12	9.208 315
2.29	3.99	0.85	Purple/black	12	9.208 316
2.54	4.24	0.85	Purple/orange	12	9.208 317
2.79	4.49	0.85	Purple/white	12	9.208 318