

Pump tubings Tygon^{®} SI 3350, with 2 colour-coded bridges

Silicone tube, platinum-cured. Meets FDA and USP Class VI and ISO 10993. 6.266 859 without bridges.

Int. Ø mm	Ext. Ø mm	Wall thickness mm	Colour	Material	PK	Cat. No.
0.64	2.44	0.90	Orange/white	Silicone	6	9.208 506
0.76	2.46	0.85	Black/black	Silicone	6	9.208 507
0.89	2.59	0.85	Orange/orange	Silicone	6	9.208 508
1.02	2.72	0.85	White/white	Silicone	6	9.208 509
1.140	2.84	0.85	Red/red	Silicone	6	9.208 510
1.3	3	0.85	Grey/grey	Silicone	6	9.208 511
1.42	3.12	0.85	Yellow/yellow	Silicone	6	9.208 512
1.52	3.22	0.85	Yellow/blue	Silicone	6	9.208 519
1.65	3.35	0.85	Blue/blue	Silicone	6	9.208 513
1.85	3.55	0.85	Green/green	Silicone	6	9.208 514
2.06	3.76	0.85	Purple/purple	Silicone	6	9.208 515
2.06	3.76	0.85	Without	Silicone	10	6.266 859
2.29	3.99	0.85	Purple/black	Silicone	6	9.208 516
2.54	4.24	0.85	Purple/orange	Silicone	6	9.208 517
2.79	4.49	0.85	Purple/white	Silicone	6	9.208 518