

Herbie tubing clamps

Beetle nylon 6.6. Simple to assemble by pressing together without the need for tools. Quick to detach with side dislocation. Good chemical resistance. Temperature resistant from -40°C to 170°C. Longterm temperature 125°C, shortterm temperature 170°C.

Min. clamping width mm	Max. clamping width mm	Clamping range	Material	PK	Cat. No.
9.2	10.3	9.2 ... 10.3 mm	Nylon	1	9.207 709
12.0	13.8	12.0 ... 13.8 mm	Nylon	1	9.207 712
13.1	15.1	13.1 ... 15.1 mm	Nylon	1	9.207 713
15.0	17.1	15.0 ... 17.1 mm	Nylon	1	9.207 715
16.9	19.1	16.9 ... 19.1 mm	Nylon	1	9.207 717
18.1	20.6	18.1 ... 20.6 mm	Nylon	1	9.207 718
20.3	23.0	20.3 ... 23.0 mm	Nylon	1	9.207 720
22.2	24.9	22.2 ... 24.9 mm	Nylon	1	9.207 722
24.5	27.6	24.5 ... 27.6 mm	Nylon	1	9.207 725

Alternativen:

Worm drive tubing clips / NORMA Distribution Center

Tubing clamps / UNEXiS

Worm drive tubing clamps / NORMA Distribution Center