

Tubing Tygon^{®} F-4040-A

Transparent, yellow, very flexible. For most fuels and industrial lubricants, fuel oils and glycol-based coolants.

- Resists swelling and cracking
- Ozone and UV light resistant

Int. Ø mm	Ext. Ø mm	Wall thickness mm	Length m	Colour	PK	Cat. No.
3.2	6.4	1.6	15	Yellow	15	6.802 073
4.8	8	1.6	15	Yellow	15	7.608 048
6.4	9.6	1.6	15	Yellow	15	6.223 959
9.6	12.7	1.6	15	Yellow	15	6.228 934
9.6	15.9	3.2	15	Yellow	15	7.972 285
15.9	22.3	3.2	15	Yellow	15	7.658 655