

# Peristaltic pump tubing E-LFL with 3 colour-coded bridges

Long flex life, high-precision tubing with 3 colour-coded bridges. Suitable for peristaltic pumps and for food and beverage dispensing. Contains no BPA or phthalates. Broad resistance to acids and alkalis. Meets USP Klasse VI and ISO 10993.

Int. Ø mm	Wall thickness mm	Length mm	Colour	PK	Cat. No.
0.89	0,85	380	Orange/orange	12	6.204 829
1.52	0,85	380	Yellow/blue	12	6.051 338
1.85	0,85	380	Green/green	12	6.056 157
2.79	0,85	380	Purple/white	12	7.661 581

## Alternativen:

Peristaltic pump tubing E-LFL with 2 colour-coded bridges / Saint Gobain

Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges / Saint Gobain