

Micropipetas monocanal Finnpiptette®; F1, volumen variable

Número de canales:

1

Capacidad	Capacidad mín. µl	Capacidad máx. µl	División µl	Precisión ± R%	Código de colores	Consejos
0,2 ... 2 µl	0,2	2,0	0,002	12,0	rosa	Finntip 10 Micro, Filter 10 Micro
0,5 ... 5 µl	0,5	5,0	0,010	6,0	rosa	Finntip 10 Micro, Filter 10 Micro
1 ... 10 µl	1,0	10,0	0,020	2,5	rosa	Finntip 10 Micro, Filter 10 Micro
1 ... 10 µl	1,0	10,0	0,020	3,5	amarillo	Finntip 250, Filter 10, Flex 20
2 ... 20 µl	2,0	20,0	0,020	3,0	turquesa	Finntip 50 Micro, Filter 50 Micro
2 ... 20 µl	2,0	20,0	0,020	3,0	amarillo	Finntip 250, Filter 20, Flex 20
5 ... 50 µl	5,0	50,0	0,200	3,0	turquesa	Finntip 50 Micro, Filter 50 Micro
5 ... 50 µl	5,0	50,0	0,200	3,0	amarillo	Finntip 250, Filter 100, Flex 20
10 ... 100 µl	10,0	100,0	0,200	3,0	amarillo	Finntip 250, Filter 100, Flex 20
20 ... 200 µl	20,0	200,0	0,200	1,8	amarillo	Finntip 250, Filter 200, Flex 20
30 ... 300 µl	30,0	300,0	1,000	1,5	naranja	Finntip 300, Filter 300, Flex 300
100 ... 1.000 µl	100,0	1000,0	1,000	1,0	azul	Finntip 1000, Filter 1000, Flex 1000
500 ... 5.000 µl	500,0	5000,0	10,000	2,0	verde	Finntip 5 ml, Filter 5 ml
1.000 ... 10.000 µl	1000,0	10000,0	20,000	2,0	rojo	Finntip 10 ml, Filter 10 ml, Flex 10 ml

Zubehör:

Soporte para micropipetas monocanal y multicanal F1 / F2 / Thermo Elect.LED