

# Micropipetas monocanal FinnpiPETTE®; F2, volumen fijo

Capacidad	Capacidad máx. µl	Precisión ± R%	Número de canales	Consejos
1 µl	1,0	4,0	1	Finntip 10, Filter 10, Flex 10, Flex Filter 10
5 µl	5,0	1,4	1	Finntip 10, Filter 10, Flex 10, Flex Filter 10
10 µl	10,0	0,9	1	Finntip 250, Filter 10, Flex 200, Flex Filter 10, 200 Ext
25 µl	25,0	0,6	1	Finntip 250, Filter 10, Flex 200, Flex Filter 10, 200 Ext
20 µl	20,0	0,6	1	Finntip 250, Filter 10, Flex 200, Flex Filter 10, 200 Ext
50 µl	50,0	0,6	1	Finntip 250, Filter 10, Flex 200, Flex Filter 10, 200 Ext
100 µl	100,0	0,4	1	Finntip 250, Filter 10, Flex 200, Flex Filter 10, 200 Ext
200 µl	200,0	0,4	1	Finntip 250, Filter 10, Flex 200, Flex Filter 10, 200 Ext
250 µl	250,0	0,4	1	Finntip 1000, Filter 1000, Flex 1000, Flex Filter 1000, 1000 Ext, Filter 1000
500 µl	500,0	0,3	1	Finntip 1000, Filter 1000, Flex 1000, Flex Filter 1000, 1000 Ext, Filter 1000
1.000 µl	1000,0	0,3	1	Finntip 5 ml, Filter 5 ml
2.000 µl	2000,0	0,3	1	Finntip 5 ml, Filter 5 ml
3.000 µl	3000,0	0,3	1	Finntip 5 ml, Filter 5 ml
5.000 µl	5000,0	0,3	1	Finntip 5 ml, Filter 5 ml
10.000 µl	10000,0	0,3	1	Finntip 10 ml, Filter 10 ml, Flex 10 ml Ext

Zubehör:  
Soporte para micropipetas monocanal y multicanal F1 / F2 / Thermo Elect.LED