

Micropipetas monocanal Finnpiquette®; F3, volumen variable

Número de canales:

1

Capacidad	Capacidad mín. µl	Capacidad máx. µl	Precisión ± R%	Código de colores	Consejos
1 ... 10 µl	1,0	10,0	2,5	rosa	Finntip 10 Micro, Filter 10 Micro, Flex 10, F
1 ... 10 µl	1,0	10,0	3,5	amarillo	Finntip 250, Filter 10, Flex 200, Flex Filter
2 ... 20 µl	2,0	20,0	3,0	turquesa	Finntip 50 Micro, Filter 50 Micro
2 ... 20 µl	2,0	20,0	3,0	amarillo	Finntip 250, Filter 20, Flex 200, Flex Filter
5 ... 50 µl	5,0	50,0	3,0	turquesa	Finntip 50 Micro, Filter 50 Micro
5 ... 50 µl	5,0	50,0	3,0	amarillo	Finntip 250, Filter 100, Flex 200, Flex Filter
500 ... 5.000 µl	500,0	5000,0	2,0	verde	Finntip 5 ml, Filter 5 ml
10 ... 100 µl	10,0	100,0	3,0	amarillo	Finntip 250, Filter 100, Flex 200, Flex Filter
20 ... 200 µl	20,0	200,0	3,0	amarillo	Finntip 250, Filter 200, Flex 200, Flex Filter
100 ... 1.000 µl	100,0	1000,0	1,5	azul	Finntip 1000, Filter 1000, Flex Filter 1000,
1.000 ... 10.000 µl	1000,0	10000,0	2,0	rojo	Finntip 10 ml, Filter 10 ml, Flex 10 ml, Flex