

Blocks for centrifuge tubes for dry block heater ECODRYBlock

Interchangeable aluminum blocks guarantee optimal thermal conductivity, uniform heating, and easy cleaning for tempering conical bottom centrifuge tubes. Designed for centrifugation of medium to large capacities for the separation of cells, particles, or biological components, as well as for the preparation and collection of pellets. Each block includes a dedicated hole to accommodate an external Pt100/Pt1000 temperature sensor, allowing direct monitoring and control of the actual sample temperature.

Dimensions (WxDxH): 95 x 76 x 51 mm

Material block: Aluminium

Type	For type of vessels	For vessels ml	No. of borings	Ø bores mm	Depth bores mm	Width mm	Depth mm	Height mm
CCT-W12-V150-A	Centrifuge tubes	15	12	17.1	44.5	95	76	51
CCT-W05-V500-A	Centrifuge tubes	50	5	29	47.6	95	76	51