

Rotors à angle fixe pour centrifugeuses Superspeed LYNX

Pour appareils:

LYNX 4000, 6000

Type	Pour r�cipients	Mat�riau du rotor	Angle �	FCR max. x-g	Vitesse max. 1/min	UC	R�f�rence
A23-6 x 100	6 x 100 ml	Aluminium	25	62976	23500	1	4.685 775
A27-6 x 50	6 x 50 ml	Aluminium	34	79057	27000	1	4.712 458
A27-8 x 50	8 x 50 ml	Aluminium	34	87207	27000	1	4.712 459
T29-8 x 50	8 x 50 ml	Titane	34	100605	29000	1	4.685 776
A22-24 x 16	24 x 16 ml	Aluminium	42	60063	22000	1	6.522 193
A21-24 x 15c	24 x 15 ml, conique	Aluminium	42	63049	21500	1	4.712 457