

# LLG-Cleaning brushes

Excellent tools for washing and cleaning of all type of laboratory bottles, flasks, tubes, condensers, burettes, etc. Brush material is manufactured from highest quality of fibre/perlon for long life cycle. Brush stem is manufactured from corrosion proof twisted steel wire. Brush handle is manufactured from polypropylene and offers extreme easy handling and use of the brush. Unique handle design enables hanging of the brush to a convenient place for easy reach. Spherical sponge top part enables access to every corner of the labware and guarantees perfect cleaning results.

Ø mm	Overall length mm	Material handle	PK	Cat. No.
12.0	235	Wire	10	6.267 298
25.0	295	Wire	10	6.267 297
40.0	380	Wire	10	6.267 296
50.0	460	Wire	10	6.267 295
70.0	460	Wire	10	6.267 294
20.0	845	Wire	10	6.267 293