

LLG-Filter crucibles, porcelain

Glazed. With porous bottom. According to DIN 12909. Temperature resistant up to 900 °C.

Max. temp.resistance:

900 °C

Material:

Porcelain

Nominal capacity ml	Ø mm	Height mm	Dimensions (ØxH)	Pore size µm	Standard	PK	Cat. No.
8	25	28	25 x 28 mm	6	DIN 12909	1	9.052 080
8	25	28	25 x 28 mm	7	DIN 12909	1	9.052 081
8	25	28	25 x 28 mm	8	DIN 12909	1	9.052 082
15	30	35	30 x 35 mm	6	DIN 12909	1	9.052 083
15	30	35	30 x 35 mm	7	DIN 12909	1	9.052 084
15	30	35	30 x 35 mm	8	DIN 12909	1	9.052 085
25	35	40	35 x 40 mm	6	DIN 12909	1	6.233 176
25	35	40	35 x 40 mm	7	DIN 12909	1	6.233 748
25	35	40	35 x 40 mm	8	DIN 12909	1	9.052 086
35	40	43	40 x 43 mm	6	DIN 12909	1	9.052 087
35	40	43	40 x 43 mm	7	DIN 12909	1	9.052 088
35	40	43	40 x 43 mm	8	DIN 12909	1	9.052 089
50	45	50	45 x 50 mm	6	DIN 12909	1	9.052 090
50	45	50	45 x 50 mm	7	DIN 12909	1	9.052 091
50	45	50	45 x 50 mm	8	DIN 12909	1	9.052 092

Zubehör:

Lids for LLG-Crucibles, porcelain / LLG Labware