

# LLG-Crimp Seals ND20, Aluminium, ready assembled

Cap size	Material	Colour	With hole	Int. bore Ø mm	Material septa	Colour septa	Thickness mm	Shore hardness	S
ND20	Aluminium	Silver	Yes	10.0	Butyl/PTFE	Grey/grey	3.00	50 A	
ND20	Aluminium	Silver	Yes	10.0	Butyl/PTFE	Grey/grey	3.00	50 A	
ND20	Aluminium	Silver	Yes	10.0	Butyl/PTFE	Grey/grey	3.00	50 A	F
ND20	Aluminium	Gold	Yes	10.0	Butyl/PTFE	Grey/grey	3.00	50 A	F
ND20	Aluminium	Silver	Yes	10.0	Silicone/PTFE	Blue transparent/transparent	3.00	45 A	
ND20	Aluminium	Silver	Yes	10.0					
ND20	Aluminium	Silver	Yes	10.0	Silicone/PTFE	White/beige	3.20	45 A	