

# LLG-Bottles, wide mouth, HDPE

Suitable for food. With carrying handle. Ideal for containing liquids, powders and granules. Very strong and durable quality for extended use over time. Screw cap can be attached to the bottle using the integrated loops. The bottles can be used by temperature from -40 &deg;C to 80 &deg;C.

Please use only in upright position.

No. of handles:

1  
-40 ... 80 &deg;C

Temperature resistance:

Nominal capacity l	Ø mm	Height mm	Dimensions (ØxH)	Neck int. Ø mm	Thread Ø mm	Min. temp. resistance °C	Max. temp. resistance °C
5	170	350	170 x 350 mm	84	94	-40	80
10	210	420	210 x 420 mm	99	109	-40	80