

# b.safe Waste Caps B 63, PE

Screw caps with free movable insert made of PTFE. Depending on the version, for connecting capillaries (O.D. 1.6/2.2/3.2 mm) and hose olives for tubing with I.D. 5.5 to 8 mm. A b.safe exhaust air filter can be screwed into the thread GL 14 (not included in the scope of supply).

## Scope of supply:

Fittings (UNF 1/4"), blind fittings (UNF 1/4"), hose olive (NPT 1/8")

Thread	Capillary connections	Hose connections	Connection for exhaust filter	No. of fittings	No. of blind plugs	No. of olives
B 63	3		1	3 x 1,6/ 2,2/ 3,2 mm	2	0
B 63	2	1	1	2 x 1,6/ 2,2/ 3,2 mm	1	1
B 63	3	1	1	3 x 1,6/ 2,2/ 3,2 mm	2	1

## Zubehör:

Fittings, PFA, for Caps and Waste Caps / Bohlender

Blind fittings, for Caps and Waste Caps / Bohlender

Exhaust Filter, PP / Bohlender

Hose connectors for Waste Caps / Bohlender