

2-way GL stopcocks, PTFE

Stopcocks with cylindrical taps and GL screwthread made of pure PTFE, giving almost total chemical resistance. Can be used at working pressures of up to 6 bar when used in combination with Bola laboratory screwthread connections (at room temperature).

A further advantage is that they can be used as expansion/reduction adapters to connect tubing of different diameters. For easy tubing connection we recommend accessory Bola laboratory screwthread connectors - please order separately.

Material: PTFE

Bore mm	Length mm	Thread	Max. temp. resistance °C	Temperature resistance	No. of nozzles	PK	Cat. No.
4	54	GL 14	110	0 ... 110 °C	2	1	9.207 153
6	64	GL 18	110	0 ... 110 °C	2	1	9.207 154
8	78	GL 25	110	0 ... 110 °C	2	1	9.207 155

Alternativen:

2-way ball valves / Bohlender