

2-way ball valves

Ball valves with GL screwthread made of chemically resistant PTFE. Maximum permissible atmospheric temperature 150°C, passing media can be as hot as 220°C.

The special construction allows the clamps to be used under pressures of up to 10 bar and also under vacuum. A further advantage is that they can be used as expansion/reduction adapters to connect tubing of different diameters. For easy tubing connection we recommend the accessory Bola laboratory screwthread connection which can be ordered separately.

Bore mm	Length mm	Thread	Material	PK	Cat. No.
3	50	GL 14	PTFE	1	9.207 305
4	80	GL 18	PTFE	1	9.207 306
8	90	GL 25	PTFE	1	9.207 307
12	100	GL 32	PTFE	1	9.207 308

Alternativen:

2-way GL stopcocks, PTFE / Bohlender