

LLG-Insertion thermometer, Type 12070, digital

Measuring range -50 to +300 °C, accuracy ± 1 °C until 150 °C, above ± 2 °C. With hold function. Probe 200 x 3.5 mm with tip made of stainless steel, with adjustable protective plastic sleeve and clip, with battery 1 x 1.5 V (LR44). Dimensions (L x W): 78 x 25 mm.

Sensor (ØxL): 3.5 x 200 mm
Battery/rechargeable battery: 1 x 1.5 V, LR44

Type	Min. temperature °C	Max. temperature °C	Temperature	Accuracy
12070	-50	300	-50 ... 300 °C	± 1 °C (<150 °C) ± 2 °C (>150 °C)

Alternativen:

- LLG-Digital pocket thermometer Type 12080 / LLG Labware
- LLG-Digital pocket thermometer Type 12090 / LLG Labware
- LLG-Insertion thermometer, Type 12060, digital / LLG Labware
- LLG-Insertion thermometer, Type 12050, digital / LLG Labware