

LLG-General purpose thermometer, red filling, immersion depth total

Graduation resistant against acids and alkalis.

Colour filling:

Red

Min. temperature °C	Max. temperature °C	Temperature	Division °C	Max. Ø mm	Length mm	Capillary	Shape therm
-35	50	-35 ... 50 °C	1	6	300	White enameled, round	Stem shape
-20	110	-20 ... 110 °C	1	6	300	White enameled, round	Stem shape
-20	150	-20 ... 150 °C	1	6	300	White enameled, round	Stem shape
-10	110	-10 ... 110 °C	1	6	260	Yellow enameled, round	Stem shape
-10	150	-10 ... 150 °C	1	6	260	Yellow enameled, round	Stem shape
-10	200	-10 ... 200 °C	1	6	300	Yellow enameled, round	Stem shape
-10	200	-10 ... 200 °C	1	6	300	White enameled, round	Stem shape
-10	250	-10 ... 250 °C	1	6	300	Yellow enameled, round	Stem shape

Alternativen:

LLG-General purpose thermometers *economy*, red filling / LLG Labware

LLG-General-purpose thermometers / LLG Labware

LLG-General-purpose thermometers, green filling / LLG Labware

LLG-General purpose thermometer, red filling / LLG Labware