

Hose tubing stopcock, HDPE

Stopcock with two-sided tubing nozzle.

- Resistant to acids and chemicals
- Pressure resistant up to 0.3 bar

Colour: Natural
Material: HDPE

Min. int. tubing Ø mm	Max. int. tubing Ø mm	Int. tubing Ø	Bore mm	Length mm	No. of nozzles	PK	Cat. No.
9	11	9 ... 11 mm	7	85	2	1	6.050 124
12	15	12 ... 15 mm	9	97	2	1	6.056 816