

Magnetic stirring bars SATELLITE 70 / 140

Stirring Bar SATELLITE 70

Extremely strong, abrasion-free and wear-free magnetic stirring bar, bearing on a tripod-construction, especially for the use in stainless steel vessels and in the case of media containing abrasive particles, e.g. sand or sediments, no direct contact between stirring bar and reaction vessel, therefore steady firm, wear-free, chemical resistance comparable to PTFE, recommended magnetic stirrers: 2mag MIX 1 XL and MIXdrive 1 XL.

Stirring Bar SATELLITE 140

Identical to SATELLITE 70, but for larger volumes/vessels, recommended magnetic stirrers: 2mag maxMIX, MAXdrive, FABdrive.

Type	PK	Cat. No.
SATELLITE 70	1	9.646 032
SATELLITE 140	1	6.263 848