

# Semi-micro balances Explorer EXR

The Explorer EXR series offers hassle-free data transfer with driverless HID (Human Interface Device). The time synchronization NTP function confirms data traceability to meet pharmaceutical regulations.

- 7-Inch touchscreen display
- Two leveling feet with digital leveling bubble
- Ventilated weighing pan
- Touchless sensors
- Full metal housing
- Removable draftshield
- Internal calibration systems
- 2 x USB-A, 1 x USB-C, 1 x Ethernet, 1 x RS232
- Optional: Bluetooth<sup>®</sup>/WiFi

Weighing plate (WxD): 80 x 80 mm  
 Dimensions (WxDxH): 230 x 393 x 370 mm  
 Weight: 7.4 kg  
 Power supply: 100 ... 240 V·AC, 50/60 Hz

Type	Max. load g	Readability mg	Reproducibility mg	Linearity	Verified calibration	Type of display	With draft shield	W m
EXR125D	82/120	0,01	0,01	&plusmn;0,06 mg		TFT	Yes	23
EXR125DM	82/120	0,01	0,01	&plusmn;0,06 mg	Verified	TFT	Yes	23
EXR125	120	0,01	0,01	&plusmn;0,06 mg		TFT	Yes	23
EXR125M	120	0,01	0,01	&plusmn;0,06 mg	Verified	TFT	Yes	23
EXR225D	120/220	0,01	0,01	&plusmn;0,06 mg		TFT	Yes	23
EXR225DM	120/220	0,01	0,01	&plusmn;0,06 mg	Verified	TFT	Yes	23
EXR225	220	0,01	0,01	&plusmn;0,06 mg		TFT	Yes	23
EXR225M	220	0,01	0,01	&plusmn;0,06 mg	Verified	TFT	Yes	23

Alternativen:  
 Semi-micro balances Explorer Plus EXP / Ohaus