

Magnetic stirring bars, PTFE, oval

For vessels with round bottom.

Ø mm	Length mm	Material sheathing	Form magnetic stirring bars	PK	Cat. No.
5.0	10.0	PTFE	Elliptical shape	10	4.688 879
6.0	15.0	PTFE	Elliptical shape	10	4.688 881
10.0	20.0	PTFE	Elliptical shape	10	4.688 882
12.0	25.0	PTFE	Elliptical shape	10	4.688 883
15.0	32.0	PTFE	Elliptical shape	10	4.688 884
15.0	34.0	PTFE	Elliptical shape	10	4.688 885
20.0	40.0	PTFE	Elliptical shape	10	4.688 886
20.0	50.0	PTFE	Elliptical shape	10	4.688 887
20.0	65.0	PTFE	Elliptical shape	10	4.688 910
20.0	70.0	PTFE	Elliptical shape	10	4.688 911

Alternativen:

Magnetic stirring bars, PTFE, octagonal / BRAND

Magnetic stirring bars, PTFE, cross shape / BRAND