

Semi-micro balances Explorer Plus EXP

The Explorer EXP series offers hassle-free data transfer with driverless HID (Human Interface Device). The time synchronization NTP function confirms data traceability to meet pharmaceutical regulations.

- 7 inch touchscreen
- Detachable terminal
- Motorised leveling system with digital leveling bubble
- Anti-pinch auto doors
- Touchless sensors
- Ionizers
- Full metal housing
- Removable draftshield
- Internal calibration systems
- 2 x USB-A, 1 x USB-B, 1 x USB-C, Ethernet, 1 x RJ11, 1 x RS232
- Optional: Bluetooth[®]/WiFi

Weighing plate (WxD): 80 x 80 mm
 Dimensions (WxDxH): 229 x 408 x 372 mm
 Weight: 7.8 kg
 Power supply: 100 ... 240 V·AC, 50/60 Hz

| Type | Max. load g | Readability mg | Reproducibility mg | Linearity | Verified calibration | Type of display | With draft shield |
|-------------|----------------|-------------------|-----------------------|-----------------|----------------------|-----------------|-------------------|
| EXP125D/AD | 82/120 | 0,01 | 0,01 | ±0,06 mg | | TFT | Yes |
| EXP125DM/AD | 82/120 | 0,01 | 0,01 | ±0,06 mg | Verified | TFT | Yes |
| EXP125/AD | 120 | 0,01 | 0,01 | ±0,06 mg | | TFT | Yes |
| EXP125M/AD | 120 | 0,01 | 0,01 | ±0,06 mg | Verified | TFT | Yes |
| EXP225D/AD | 120/220 | 0,01 | 0,01 | ±0,06 mg | | TFT | Yes |
| EXP225DM/AD | 120/220 | 0,01 | 0,01 | ±0,06 mg | Verified | TFT | Yes |
| EXP225/AD | 220 | 0,01 | 0,01 | ±0,06 mg | | TFT | Yes |
| EXP225M/AD | 220 | 0,01 | 0,01 | ±0,06 mg | Verified | TFT | Yes |