

# Blender bags BagSystem<sup>&reg;</sup>, instaBAG<sup>&reg;</sup> Fraser &frac12; with dehydrated medium

With medium for primary selective enrichment or as a diluent for the detection and enumeration of *Listeria monocytogenes* and *Listeria* spp. (ISO 11290-1, ISO 11290-2) in food and environmental samples.

Nominal capacity ml	Width mm	Height mm	Dimensions (WxH)	Colour	PK	Cat. No.
225	190	300	190 x 300 mm	Transparent	80	6.316 572