

# Staining jar, PMP (TPX<sup>&reg;</sup>), Schiefferdecker pattern

Accepts 20 slides 76mm x 26mm (10 pairs back to back). Transparent. Complete with lid.

Note: Xylene can cause etching after 2 days at room temperature.

With lid:

Yes

Width mm	Depth mm	Height mm	Dimensions (WxDxH)	Colour	Material	PK	Cat. No.
86	70	51	86 x 70 x 51 mm	Clear	PMP	1	9.161 352

Zubehör:

Staining jar, glass, Coplin / Glaswarenfabrik Karl Hecht

Slide holder, ABS / Kartell S.p.A.

Alternativen:

Staining jar, PMP (TPX<sup>®</sup>), Hellendahl pattern / Kartell S.p.A.

Staining troughs / Glaswarenfabrik Karl Hecht

Ethidium Bromide Extractor / GVS Microfiltrazione SRL

Staining jar, PMP (TPX<sup>®</sup>) / Kartell S.p.A.