

# Spectrophotometer WPA Biowave 3 Wireless, with Bluetooth

Scattered light: <0.5 %T (340 nm)  
Illuminants: Xenon  
Interfaces: USB, Bluetooth  
Dimensions (WxDxH): 260 x 390 x 100 mm  
Weight: 3 kg  
Power supply: 100 ... 240 V·AC, 50/60 Hz  
Warranty: 2 years

Type	Wavelength nm	Accuracy	Spectral bandwidth nm	Photometric measuring range	Measuring method	Width mm
Biowave 3 Wireless	190 ... 1100	±2 nm	5	-0.300 ... 2.500 A0.3 ... 199 %T	UV-VIS	260

## Alternativen:

Spectrophotometer WPA Biowave 3+ Wireless, with Bluetooth / Biochrom Ltd.  
Spectrophotometer WPA Lightwave 3 Wireless, with Bluetooth / Biochrom Ltd.  
Spectrophotometer WPA Lightwave 3+ Wireless, with Bluetooth / Biochrom Ltd.  
Spectrophotometer WPA Biowave 3+ CT Wireless, with Bluetooth / Biochrom Ltd.  
Spectrophotometer WPA Lightwave 3+ CT Wireless, with Bluetooth / Biochrom Ltd.