

Membrane filters, Nylon

Hydrophilic, durable and compatible with a wide range of aqueous solutions, organic solvents and for the filtration of biological preparations and media. Autoclavable at 121 °C.

Ø mm	Pore size µm	Material	Autoclavable	PK	Cat. No.
13	0.20	Nylon	Yes	100	9.951 042
13	0.45	Nylon	Yes	100	9.951 045
25	0.20	Nylon	Yes	100	6.240 949
25	0.45	Nylon	Yes	100	9.951 046
47	0.20	Nylon	Yes	100	9.951 043
47	0.80	Nylon	Yes	100	9.951 049
47	0.45	Nylon	Yes	100	9.951 047
90	0.20	Nylon	Yes	50	9.951 044
90	0.45	Nylon	Yes	50	9.951 048