

Aspirator Carboys behroplast[®], HDPE

With welded stopcock. Food-safe. Leakproof, even after prolonged use.

With stopcock: Yes

Nominal capacity l	Width mm	Depth mm	Height mm	Dimensions (WxDxH)	Material	Food safe	PK	Cat. No.
5	160	190	240	160 x 190 x 240 mm	HDPE	Yes	1	9.140 100
10	190	190	360	190 x 190 x 360 mm	HDPE	Yes	1	9.140 101
20	250	280	370	250 x 280 x 370 mm	HDPE	Yes	1	9.140 098
25	240	300	450	240 x 300 x 450 mm	HDPE	Yes	1	9.140 099
30	260	300	450	260 x 300 x 450 mm	HDPE	Yes	1	9.140 102
60	330	390	670	330 x 390 x 670 mm	HDPE	Yes	1	9.140 103